Document No.: 472313373
Revision: A
Change Order: CO-035068
Title: PRODUCT INSERT PLASMACELL-C DISPOSABLE SET FOR USE WITH SPIKESMART SYSTEM K6R2256
Page: 1 of 4
Last Modified: 04 Jan 2019 JML
Color: Black

PROOFREADING INSPECTION / RELEASED ARTWORK

Proofreading Approval:

Print Name _____

Signature _____

Date _____

Code REF K6R2256

60x

(en) PLASMACELL-C Disposable Set for Use with SPIKESMART System

For Collection of Virtually Cell Free Plasma by Membrane Filtration.
For Use with the AURORA Instrument or the AUTOPHERESIS-C Instrument.

Store at Room Temperature. Protect from freezing. Avoid excessive heat.

Definition of "Room Temperature":
"The temperature prevailing in a working area."

Reference: United States Pharmacopeia, General Notices.
United States Pharmacopeial Convention, Inc.
12601 Twinbrook Parkway, Rockville, MD.

Description:

The PLASMACELL-C disposable set is a single unit open system consisting of the PLASMACELL-C separation device (containing the membrane filter), the reinfusion reservoir, two (2) transducer protectors, two (2) solution spikes, one (1) female Luer for plasma container connection, and one (1) male Luer for needle connection.

The extracorporeal volume of the set, including the separation device and the reinfusion reservoir, is approximately 200 mL.

Indications and Usage:

The system is an automated donor plasmapheresis system designed to collect plasma by membrane filtration. Plasma must be handled (freezing, storage conditions, labeling, etc.) per center requirements and applicable local regulations.

Contraindications:

The use of the AURORA and AUTOPHERESIS-C plasmapheresis systems are contraindicated in those cases where adequate anticoagulation cannot be achieved.

Precautions:

1. Carefully follow the installation instructions provided in the Operator's Manual to avoid misconnection and leaks.
2. Before operating the AURORA or AUTOPHERESIS-C system, the operator should be trained and knowledgeable in the system components, controls, and operating protocol. Although the system is automated, operator surveillance is important during plasmapheresis.
3. The suitability of the donor for this procedure and frequency of donation should be determined according to center requirements and the requirements of the U.S. Food and Drug Administration, or according to local regulatory requirements, as appropriate outside the U.S.
4. For any deviation in performance, refer to the TROUBLESHOOTING section of the Operator's Manual.

As with any procedure involving extracorporeal circulation, there is a risk of hemolysis, air embolism or blood clotting. To avoid these risks, do not deviate from the AURORA Operator's Manual or the AUTOPHERESIS-C Operator's Manual. In the event that color change is observed in the plasma line, follow the directions in the TROUBLESHOOTING section of the Operator's Manual.

During connection of solutions, do not interchange the anticoagulant and saline lines.

Adverse Effects:

Adverse donor reactions identical to those seen with routine blood collection procedures may occur, which include:

- Lightheadedness
- Fainting
- Vomiting
- Hyperventilation
- Hematoma formation at the site of phlebotomy

Reactions unique to plasmapheresis collection procedures may also occur. Reinfusion of saline or donor blood may induce chills if the procedure is carried out at very cool temperatures. Infusion of anticoagulants containing citrate at rapid rates may lead to donor symptoms of modest hypocalcemia due to chelation of calcium by unmetabolized citrate. Such reactions are usually manifested by a tingling feeling, often around the mouth of the donor. Return of the donor blood should be temporarily slowed or halted if such symptoms occur.

Symptoms of severe hypocalcemia, although rare, include tetany, convulsions, cardiac arrhythmia, and death.

Improper operating conditions may result in complications such as blood loss, hemolysis, air embolism, and blood clotting.

This container is of PVC material that used di-(2-ethylhexyl) phthalate (DEHP) as plasticizer, where DEHP is reported to affect the development of male sexual organs and sperm formation in young animal tests. Although infinitesimal amount of DEHP may be eluted in such PVC containers, there is very little or no risk from exposure to DEHP. Hence, as the risk that may appear from not using this medical device is much larger than the risk concerned due to DEHP, there is no need to avoid using it.

Document No.: **472313373**
 Revision: **A**
 Change Order: **CO-035068**
 Title: **PRODUCT INSERT PLASMACELL-C DISPOSABLE SET FOR USE WITH SPIKESMART SYSTEM K6R2256**

Page: 2 of 4
 Last Modified: 04 Jan 2019 JML
 Color: Black

PROOFREADING INSPECTION / RELEASED ARTWORK		
Proofreading Approval:		
Print Name _____	Signature _____	Date _____

Directions for Use:

1. Install the separation device and the reservoir.
 - A. Separation Device Installation: Insert the separation device in the separation device motor cup, slowly raise device support, rotate separation device until it has clicked and locked on to device support.
 - B. Reinfusion Reservoir Installation: Insert the reservoir ridge into the channel, push upward and then press bottom to locks the reservoir into the channel. This channel has four sensors which measure the height of blood hence reservoir should be mounted correctly.
2. Place plasma line in the Hb detector and close door. Position the plasma line into the plasma clamp and tube guide. Aseptically connect a plasma collection container to the plasma line and hang the container on the weigh scale. (collection container must hang unobstructed)
3. Thread cell line through cell pump by placing tubing into tube guide at the top and bottom of the pump.
4. Position saline line into the saline clamp.
5. Position reinfusion line into reinfusion clamp and blood line into blood clamp. Load blood line into blood pump by placing tubing into tube guide at the top and bottom of the pump. Then insert the donor line into the air detector and turn the wheel to close.
6. Install both pressure lines to the pressure transducer assembly: Open the transducer cover, install both lines by looping it behind the pressure transducer cover, pull lines firmly until tubing can be felt sliding into the pressure transducer cover slots. Connect lines to their respective sensor ports and push cover closed.
7. Load anticoagulant line into anticoagulant pump by placing tubing into tube guide at the top and bottom of the pump.
8. The lines installed into the pumps should be centered on pump rollers without any slack after the pump cover is closed. After installing the set, recheck all connections are secure.

Refer to the Operator’s Manual for the AURORA or AUTOPHERESIS-C system for the installation procedure and complete instructions for use.

Dispose of waste in appropriate biohazard container or according to local regulatory requirements.

Symbols with Definitions:

- | | | | |
|---|---|---|--|
|  | Caution, consult instructions for use |  | Lot |
|  | Sterilized by irradiation. Sterile fluid path. |  | Expiration date |
|  | Non-pyrogenic fluid path |  | Manufacturer |
|  | Do not reuse |  | Manufacturing facility/
Manufactured by |
|  | Do not use if the product sterile barrier system is compromised |  | Fragile |
|  | Contains or presence of phthalate: Diethylhexylphthalate (DEHP) |  | This way up |
|  | Code |  | Recyclable |

Fresenius Kabi AG
 61346 Bad Homburg / Germany
 Tel.: +49 (0) 61 72 / 686-0
 www.fresenius-kabi.com

Fenwal International, Inc.
 Carretera Sanchez Km 18.5
 Parque Industrial Itabo
 Zona Franca Ind. De S.C.
 Haina, Dominican Republic

© 2019 Fresenius Kabi AG. All rights reserved.

Document No.: 472313373
Revision: A
Change Order: CO-035068
Title: PRODUCT INSERT PLASMACELL-C DISPOSABLE SET FOR USE WITH SPIKESMART SYSTEM K6R2256
Page: 3 of 4
Last Modified: 04 Jan 2019 JML
Color: Black

PROOFREADING INSPECTION / RELEASED ARTWORK

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