



## Agilia SP TIVA Agilia SP TIVA WiFi

Syringe infusion pumps providing  
intuitive intravenous anaesthesia

Agilia, Vigilant and Centerium are registered trademarks by Fresenius Kabi in selected countries. Due to our policy of continuous product development as well as changes in standards, the features described are subject to change. Please contact us for the most updated information. These medical devices are regulated by health authorities. Please refer to the respective CE Declaration of Conformity for the applicable Medical Device Directive or Regulation. Pictures are not contractual.

Dose rate programming for lower risk of errors in prescriptions

TCl infusion mode

Wide flow rate ranging from 0.1-1200 mL/h and an accuracy of +/- 3 %

Adaptable to all protocols

Up to 19 drug libraries

Easy and intuitive to use

Stand alone or rack mounted

2 pressure modes and Dynamic Pressure System (DPS)

# Agilia SP TIVA & Agilia SP TIVA WiFi

## Infusion

### Flow rate range

0.1-1200 mL/h  
0.1 mL/h increment from 0.1 to 99.9 mL/h (optionally 0.01 mL/h from 0.10 to 9.99 mL/h), 1 mL/h from 100 to 1200 mL/h  
Flow rate can be restricted for each drug (soft and hard limits) with Vigilant Master Med, The Dose Error Reduction Software

### Flow rate accuracy

+/- 1% on mechanism ; +/- 2% with syringes

### Syringes capacities

5, 10, 20, 30/35, 50/60 CC

### Type of syringe

Up to 100 types

### Infusion Modes

- mL/h mode
- Dose rate modes: ng/h, ng/kg/min, ng/kg/h, microg/min, microg/h, microg/kg/min, microg/kg/h, mg/min, mg/h, mg/24h, mg/kg/min, mg/kg/h, mg/kg/24h, mg/m<sup>2</sup>/h, mg/m<sup>2</sup>/24h, g/h, g/kg/min, g/kg/h, g/kg/24h, mmol/h, mmol/kg/h, mmol/kg/24h, mU/min, mU/kg/min, mU/kg/h, U/min, U/h, U/kg/min, U/kg/h, kcal/h, kcal/24h, kcal/kg/h, mEq/min, mEq/h, mEq/kg/min, mEq/kg/h
- TCI infusion mode
- Dilution setting: -- units / mL or -- units / -- mL
- With or without loading dose
- Volume or dose / time: 0.1 - 99.9 mL; 00 h 01 - 96 h 00
- Volume limit: 0.1 - 999 mL

### TCI Mode

- Pharmacokinetic Models:
- Marsh & Schnider with Plasma and TCI Effect target for Propofol administration in adults
  - Kataria & Paedfusor with Plasma target for Propofol administration in pediatrics
  - Minto, Gepts & Scott for with Plasma and TCI Effect target for Remifentanil, Sufentanil and Alfentanil administration in adults

### Volume/Dose infused

Volume: 0.1 - 999 mL / Dose: 0.1 - 9999 units

### Priming

3 modes: mandatory, not mandatory, or advised / Rate: 1200 mL/h

### Bolus

- Direct bolus: Rate: 50 - 1200 mL/h (50 mL/h increment).
- Programmed bolus (dose or volume / time): 0.1 - 99.9 mL 0.01 - 9999 units / 1 second - 24 h

### Induction dose

Dose / time: 0.01 - 9999 units / 1 second - 24 h  
Rate auto-calculation

### End infusion (V/T & VL)

KVO: adjustable from 0.1 to 5 mL/h, continuous infusion or stop

### Quick start

Not mandatory prime set by default resulting in fast start if user does not prime with bolus button ensuring programmed flow rate to be reached faster

### Pause

Programming from 1 minute to 24 hours, increments from minute to minute

### Data log event

1500 data log events in real time

### Graphical history

Volume / dose infused, pressure, flow rate

### Night mode

The night mode decreases the brightness of the screen and the green lights. The key press sound can optionally be turned off. The night mode can be programmed manually or automatically in a variable time range

### Profiles

Basic profile infusion without any display of the drug names and with TCI modes.  
19 Custom profiles configurable with Vigilant Master Med: Drug library to be created with up to 200 drugs  
Drug list with a list of drug names  
Configuration only: custom pump configuration without Drug Name

## Pressure management

### Pressure modes

2 modes available: variable or 3 pre-set levels - Range from 50 to 900 mmHg. (25 mmHg increment from 50 to 250 mmHg / 50 mmHg increment from 250 to 900 mmHg). Can be enabled / disabled and adjusted

### DPS

The Dynamic Pressure System - DPS - warns of pressure variations. A risk of obstruction or a possible leak in the infusion line can thus be anticipated

### Pressure monitoring

Graphic representation of the pressure in the infusion line and of the pressure limit thanks to the pictogram

### Anti-bolus system

Reduces significantly bolus after occlusion release\*  
≤0.35 mL max for a 50 mL syringe  
\* Test condition: back pressure: 0 mmHg - Syringe: BD Precise

## Alarms / Pre-alarms / Security

### Pump status

GREEN for infusion in progress, ORANGE for Low and Medium priority, RED for High priority - Visibility at 4 meters minimum  
All alarms are expressed by means of light indicators, written words, pictograms and sound beeps

### Syringes installation control

Syringe barrel clasp check, plunger head detection, anti-siphon system check, flange detection

### Infusion control

Occlusion pressure pre-alarm, occlusion pressure alarm, end of infusion pre-alarm, end of infusion alarm, volume limit pre-alarm, volume limit alarm, hard and soft flow rate limits, start infusion at pause end

### Device control

Disengaged driving mechanism alarm, low battery pre-alarm, discharged battery alarm, battery capacity display in hours and minutes, unconfirmed programming, technical malfunction alarm (auto-test, rotation), drive system advance check, watchdog check, communication connection failure, plug-head disengagement, auto-lock/lock code (on Keypad)

### Maintenance

Preventive maintenance warning

## Technical specifications

### Manual pusher

Protection for the ongoing infusion thanks to "Push-Guard"

### Display

Blue graphic LCD monochrome, size 66 mm x 33 mm (256 x 128 pixels)

### Swinglock clamp

Versatile clamp that allows the fixation on a rail or on a pole (Pole: 20-40 mm max. / Rail: 25-35 x 10 mm)

### Stackability

Up to 3 devices self-stackable on a pole

### Dimensions (H/W/D) / Weight

135 x 345 x 170 mm / Approximately 2.1 kg

### Battery

Characteristics: 7.2 V 2.2 Ah - Li-ion Smart battery, remaining battery life and battery charge level available on the display

Battery Life (when fully charged):

- Agilia SP TIVA and Agilia SP TIVA WiFi (WiFi disabled / not used):

> 11 h at 5 mL/h

- Agilia SP TIVA WiFi (WiFi enabled):

> 6 h at 5 mL/h

Battery recharge:

- Pump OFF: < 6 h

- Pump ON: < 20 h

### Waterproofness

IP22

### Power supply

100 V - 240 V AC - 50 / 60 Hz with functional earth

## Compliance

### Electromagnetic compatibility EMC

IEC 60601-1-2, IEC 60601-2-24

### CE Mark

Complies with the 93/42/EEC Medical Device Directive or the Medical Device Regulation (EU) 2017/745

### Electrical Compliance

Protection against leakage current: defibrillation-proof type CF  
Protection against electric shocks: class II in accordance with IEC 60601-1

### Alarm system

IEC 60601-1-8

### Usability Engineering

IEC 60601-1-6 and IEC 62366

## Wireless LAN (For Agilia SP TIVA WiFi only)

### Technology

IEEE 802.11 a/b/g/n. Frequency Band: 2.400 -> 2.500 GHz (2.4 GHz is ISM band) / 4.900 -> 5.850 GHz (High Band)

### Modulation

OFDM with BPSK, QPSK, 16-QAM, and 64-QAM 802.11b with CCK and DSSS

### Wireless Security

WPA/WPA2-Enterprise, WPA/WPA2-PSK

### Network Protocols

TCP, IPv4, DHCP, HTTP  
CE, FCC and IC compliant