


#### How does Prepare Plus® PN work?

Just enter the nutritional needs of the patient in terms of either the products or the components. Then click "check stability" and Prepare *Plus*® PN will indicate if your PN regimen is compatible and for how long.



Should your regimen be unstable, Prepare *Plus*® PN will help you to customise your formulation to obtain stability. Prepare *Plus*® PN offers suitable advice and furthermore has an intuitive interface with a red-green scale indicating your current position within the tested range to simply edit the admixture!



Now you can rest assured - you can rely on admixtures found to be compatible and stable by Prepare Plus® PN!



## What do I need to be prepared for Prepare $Plus^{\circ}$ PN?

You just need a computer running with Windows 2000 SP4 or higher, an internet connection/USB port for updates and a printer. The program includes a help functionality to guide you along the way.






Fresenius Kabi Deutschland GmbH 61346 Bad Homburg v.d.H. Germany

Phone: +49 61 72 / 686-0 Fax: +49 61 72 / 686-7253 www.fresenius-kabi.com

# PARENTERAL NUTRITION



## Compatibility & Stability Software











Admixtures for parenteral nutrition (PN) are very complex in terms of physical and chemical compatibility as they contain a broad range of components in various amounts and numerous combinations. Assessing this is complex and requires time, knowledge and experience. Fresenius Kabi takes your request for reliable and extensive compatibility data


seriously and continuously investigates new compositions. To make this data available to you, Prepare *Plus*® PN has been developed. At your finger tips, you will have a high performance, ready to use, intelligent software, setting new standards for PN compounding.




#### What are the key features of Prepare Plus® PN?

Prepare Plus® PN makes admixing easier by:

- Consistently checking for compatible and stable admixtures
- Confirming container compatibility
- Offering specific advice to amend unstable admixtures to make them stable
- Providing stability statements, work sheets, and even labels







- Available whenever you need it
- Immediate access to updates
- Easy to use with an intuitive user-friendly interface
- Supports both adult and paediatric PN regimens



## What makes Prepare Plus® PN so exceptional?

- GMP-validated and certified software
- Extensive database compiled over many years
- All results are based on actual laboratory stability tests\*
- Intelligent guidance based on a unique search algorithm





#### How will Prepare Plus PN help you to save time?

- Immediate access to reliable compatibility data
- Clear advice on the process of mixing
- Saving and recall of standard compositions
- Display of final composition per kg or per litre
- Automatic production of stability statements, work sheets and labels
- Support tools for concentration calculations, conversion of units and volume overages



### Where did the data in Prepare Plus® PN come from?

Rather than extrapolations or stability equations, all stable results given by Prepare *Plus*® PN are based on actual laboratory stability tests using the well-established Fresenius Kabi range concept. This approach gives maximum flexibility for the addition of all admixture components.

\*Data on file