



**FRESENIUS
KABI**

caring for life



CompoSeal[®] Mobilea II Portable Tube Sealer

The CompoSeal Mobilea II is a lightweight, portable radio frequency (RF) sealer designed to seal medical PVC tubing in both mobile donation settings and processing environments.

The power management system uses high capacity Lithium Polymer battery technology to combine quality sealing performance and fast charging.

High capacity Lithium Polymer battery

Portable, lightweight design

Built for demanding use in mobile settings

Consistent seal quality

Disclaimer: This document is intended for an international audience. Please contact your local sales representative for availability in Canada. For questions or product concerns please contact Fresenius Kabi Canada (165 Galaxy Blvd, Suite 100, Toronto, ON, M9W 0C8) at 1-877-779-7760 or visit the website at www.fresenius-kabi.ca.

CompoSeal Mobilea II

Portable Tube Sealer

Portable lightweight sealer with high battery capacity

- Lithium Polymer battery provides up to 1,000 seals when fully charged (validated using common sized blood bag tubing)
- Designed for portable use
 - Lightweight 1.8 lb power pack
 - Robust aluminum housing
 - Protective padding for high-demand environments
- Recharge an empty battery within 2 hours
- Consistent, quality seals
- Quick and easy operation
- Automatically adjusts sealing time based on tube thickness



Description	Weight	Dimensions (W x D x H)
CompoSeal Mobilea II set	5.36 lbs	15.4 x 4.3 x 11.8 in
Includes:	1.8 lbs	4.3 x 2.1 x 6.9 in
Powerpack (battery and generator)	2.6 lbs (weight in use)	8.2 x 3.5 x 1.3 in
Hand-held sealer		0.7 lbs
Cable	0.1 lbs	57.4 x 0.3 x 0.3 in
Charger	0.5 lbs	2.5 x 1.5 x 4.1 in
Transport case	2.0 lbs	15.3 x 4.3 x 11.8 in

For more information please contact ProductSupport@fresenius-kabi.com or 800-333-6925.

Please refer to Operating Instructions for a complete list of warnings and precautions associated with the use of this device.

CompoSeal is a registered trademark of Fresenius Kabi AG.
Copyright © 2019 Fresenius Kabi AG. All rights reserved.