



**FRESENIUS
KABI**

caring for life



CompoCool[®] WB

Cooling System

Improved cooling system for standardized, controlled, secured storage and transport of donated whole blood.

Insulated transport boxes

Small cooling element with integrated temperature control

Ergonomic box (washable), cartridge and trolley

Secure blood transport

CompoCool[®] WB

Improved cooling system for safe and standardized cooling of whole blood at 20 - 24°C

CompoCool WB's closed and lockable transport box with integrated insulation against warm and cold environmental conditions.

- Allows a controlled cooling process and continuous cooling chain also for long term storage
- Enables cooling at 22°C (± 2°C) within 2 hours for up to 24 hours
- Small cooling element with temperature and shelf life indication for fail-safe operation

Improved filtration

- Improves quality and blood processing of leukocyte depletion, increased and reliable PTL-yield, stabilization of factor VIII during storage (Vox. Sang. 1989; 56:145-150)

Easy to clean

- Suitable for washing in a dishwasher
- Removable lid

Optimal for mobile use

- Small, stackable box for 6 blood bags enables ergonomic carrying and easy handling
- Foldable trolley, easy to carry



Each cooling element is equipped with a serial number for shelf life indication and barcode. Clear and easy instruction with pictures, attached to the box lid.



Small and light cooling element optimizes handling and storage in the refrigerator (min. 16 h at 4°C). A thermo foil on each cooling indicates the actual temperature.



Transportation box capacity for 6 filled in-line blood bag systems.

The new cooling system for protected whole blood during transport and storage

Ergonomic handling and easy to move and store, the trolley requires minimum space with safe transportation of up to four full boxes



Sealing with plastic plugs (seals) warrants the cooling chain according to the validated temperatures.

CompoCool® WB Cooling System

Ordering information

For more information such as literature, technical details, working instructions and equipment, please contact your local sales representative.

Process efficiency to help you achieve more

CE marked, Class IIa according Directive 93/42/EEC as amended 2007/47/EC



**FRESENIUS
KABI**

caring for life

Fresenius Kabi AG
61346 Bad Homburg
Germany

Phone: +49 61 72 608-0
Fax: +49 6172 608-5786

compocool@fresenius-kabi.com
www.fresenius-kabi.com

