



**FRESENIUS
KABI**

caring for life



CompoGuard

Whole Blood Collection System

Optimize your whole blood collections

The complete line of CompoGuard Whole Blood Collection Systems is designed to help optimize whole blood collections in multiple collection environments.

In-line tray for optimal bag positioning and accurate weighing

Visual and audible donation prompts for improved blood collection

Safety clamp prevents air from entering collection containers

Compact size and durable transport case make CompoGuard ideal for mobile donation settings

Disclaimer: This document is intended for an international audience. Please contact your local sales representative for availability in Canada. For questions or product concerns please contact Fresenius Kabi Canada (165 Galaxy Blvd, Suite 100, Toronto, ON, M9W 0C8) at 1-877-779-7760 or visit the website at www.fresenius-kabi.ca.

CompoGuard

Whole Blood Collection System

The complete line of CompoGuard Whole Blood Collection systems is designed to help optimize whole blood collections in multiple collection environments.



	CompoGuard Basic	CompoGuard Data	CompoGuard Complete
Description	<p>CompoGuard Basic includes key functions, including accurate weighing and mixing, adjustable collection volumes and low flow alarms, for blood collections at an affordable cost.</p> <ul style="list-style-type: none"> • In-line tray • Integrated safety clamp • Large display • Easy-to-clean soft keys 	<p>CompoGuard Data includes all the functionality and portability of CompoGuard Basic, plus:</p> <ul style="list-style-type: none"> • DECT* option and USB port for collection data transfer • Optional: DonationMaster Net Data Management software system <p>*DECT=Digitally enhanced communication technology</p>	<p>CompoGuard Complete offers all the features of the CompoGuard Data, plus:</p> <ul style="list-style-type: none"> • Sealer and rechargeable battery • Gooseneck designed to improve work station for operator comfort • Basic Model + DECT* option + Gooseneck + USB port + sealer + rechargeable battery <p>(Model shown with optional scanner)</p>
Dimensions (H x W x D)	6.5 X 8.8 X 17.5 in	6.5 X 8.8 X 17.5 in	36 X 8.8 X 17.5 in Height is with Optenna extended
Weight	7.5 lbs	10.5 lbs	18 lbs



Optional Transport Case

Dimensions (H x W x D): 11 X 9.8 X 20.8 in

Weight: 6.17 lbs

Charge CompoGuard while stored, using the easy-access charging port

For more information please contact ProductSupport@fresenius-kabi.com or 800-333-6925.

For complete instructions for use and a list of warnings and precautions, refer to the Operating Instructions.

Copyright © 2019 Fresenius Kabi AG. All rights reserved.