



Amika® Enteral Sets


Overview of ENFit and Universal Enteral Pump Sets




Copyright(C) 2024 Fresenius Kabi AG. Reproduction and distribution - only with prior written authorization

7-393098 (06.18/AC)

Overview of ENFit and Universal Enteral Pump Sets



Amika⁺



GraviSet Bag, ENFit (7751938)



Enteral Giving Set Extension (7755693)

	EasyBag		VarioLine		Bag		Bag Mobile	
ART. NR.	7751907	7751900	7751909	7751919*	7751715T	7751914*	7751913	7751906
UPPER END								
ENPlus spike	●	●	●	●				
VarioLine bottle adaptor			●	●				
Attached Bag					●	●	●	●
LINE COMPONENTS								
Drip chamber	●	●	●	●	●	●		
Roller clamp							●	●
ENFit medication port	●	●	●	●	●	●	●	●
Cover on medication port				●		●		
LOWER END								
Universal Connectors (adaptor for ENLock and Luerlock)	●		●		●		●	
ENFit connector	●	●	●		●		●	●
ENFit connector with cover				●		●		
Tubing length (cm)	250	250	250	250	250	250	242	242
Length med. port to lower end (cm)	40	40	40	40	40	40	40	40

Not made with DEHP

* Pump set with transparent cover on the lower-end ENFit connector and tethered cover on the ENFit medication port